
import numpy as np import pandas as pd import tensorflow as tf import gdcm import pylibjpeg import scipy.ndimage import matplotlib.pyplot as plt import matplotlib.patches as patches import os import nibabel as nib from glob import glob import sys import cv2 as cv

The goal is to identify fractures in CT scans in the cervical spine at the same level of accuracy as radiologists.



Specifically, we are looking at C1-C7. We also want the overall probability any fractures in the cervical spine.

The submission will be evaluated using a weighted multi-label logarithmic loss. The binary weighted log loss function for label j on exam i is specified as:

$$L_{ij} = -w_j \left(y_{ij} \log(p_{ij}) + (1-y_{ij}) \log(1-p_{ij})
ight)$$

There are also some weights associated with this:

Category: Weight

Vertebrae negative: 1 Vertebrae positive: 2 Patient negative: 7 Patient positive: 14

Competition Link

```
In [7]: train = pd.read_csv("../rsna-2022-cervical-spine-fracture-detection/train.csv")
    train_bb = pd.read_csv("../rsna-2022-cervical-spine-fracture-detection/train_bounding_boxes.csv")
    id_num = "1.2.826.0.1.3680043.1363"

    train.loc[train["StudyInstanceUID"] == id_num]
#Some of the rows have an indication where the spine is fractured
#Select a specific patient ID to look at the bounding boxes + where the fracture is
```

Out[7]: StudyInstanceUID patient_overall C1 C2 C3 C4 C5 C6 C7

4 1.2.826.0.1.3680043.1363 1 0 0 0 0 1 0 0 0 0

A dicom file consists of a header (demographic info for patient + acquisition, dimensions, matrix size) and a single attribute for pixel intensity data

Preamble (128 bytes)

Prefix - 'D', 'I', 'C', 'M'

Header:

Data Set

In [8]: train["StudyInstanceUID"].nunique()

- Group 1 (0002)
 - Element 1 (0002,0000)
 - Element 2 (0002,0001)
 - Element 3...etc.
- Group 2 (0008)
- Group 3...etc.

Image Pixel Intensity Data:

We know that these files are ≤ 1 mm slice thickness. We also have axial orientation and bone kernel.

```
#Number of unique study instances (aka patient images)

Out[8]: 2019

In [7]: #Looking at an example of the meta data included img = pydicom.dcmread(f"../rsna-2022-cervical-spine-fracture-detection/train_images/{id_num}/120.dcm") print(img)
```

```
(0002, 0001) File Meta Information Version OB: b'\x00\x01' (0002, 0002) Media Storage SOP Class UID UI: CT Image Storage
  (0002, 0003) Media Storage SOP Instance UID UI: 1.2.826.0.1.3680043.1363.1.120
 (0002, 0010) Transfer Syntax UID

(0002, 0012) Implementation Class UID

(0002, 0013) Implementation Version Name

UI: Explicit VR Little Endian
UI: 1.2.40.0.13.1.1.1

UI: 1.2.40.0.13.1.1.1
   -----
(0008, 0023) Content Date DA: '20220727'
(0008, 0033) Content Time TM: '183912.555083'
(0010, 0010) Patient's Name PN: '1363'
(0018, 0050) Slice Thickness DS: '1.0'
(0020, 0004) Study Instance UID UI: 1.2.826.0.1.3680043.1363
(0020, 000e) Series Instance UID UI: 1.2.826.0.1.3680043.1363
(0020, 00032) Image Position (Patient) DS: [-149.2080078125, -350.2080078125, 173]
(0020, 0037) Image Orientation (Patient) DS: [1, 0, 0, 0, 1, 0]
(0028, 0004) Photometric Interpretation CS: 'MONOCHROME2'
(0028, 0011) Columns
(0028, 0030) Pivel Secrica
  (0008, 0018) SOP Instance UID UI: 1.2.826.0.1.3680043.1363.1.120
  (0028, 0030) Pixel Spacing
                                                                        DS: [0.583984375, 0.583984375]
  (0028, 0100) Bits Allocated
                                                                        US: 16
  (0028, 0101) Bits Stored
                                                                         US: 12
US: 11
  (0028, 0102) High Bit
  (0028, 0103) Pixel Representation
                                                                    US: 0
  (0028, 1050) Window Center
                                                                        DS: [450, 40]
  (0028, 1051) Window Width
                                                                        DS: [1500, 350]
                                                                      DS: '-1024.0'
  (0028, 1052) Rescale Intercept
                                                                           DS: '1.0'
  (0028, 1053) Rescale Slope
  (7fe0, 0010) Pixel Data
                                                                           OW: Array of 524288 elements
```

We can now see the file meta information derived from the header of this dicom image. The tag "Media Storage SOP Instance UID" (in this case: UI: 1.2.826.0.1.3680043.14.1.1) and "Transfer Syntax UID" tell us that this image can unpacked using lossless JPEG conversion. We need GDCM and pylibjpeg

Image Position (Patient) DS: [-73.87109375, -222.37109375, 620]


Dataset.file_meta -----

Gives the x-axis, y-axis, and the z-axis. Z-axis shows the position in the sagittal plane, but the image is shown in the axial plane

We also know that the images are 512 x 512, that will impact the preprocessing for modelling. We need to calculate the voxel size + transform the data into Hounsfield Units since a CT-scan tells us about the radiodensity of an object

```
In [8]: bound_box = (train_bb[train_bb["StudyInstanceUID"] == id_num]).reset_index()
In [15]: #https://github.com/pydicom/pydicom/blob/master/examples/image_processing/reslice.py
         paths = glob(f"../rsna-2022-cervical-spine-fracture-detection/train_images/{id_num}/*.dcm")
         files = [pydicom.dcmread(path) for path in paths]
         images = [pydicom.pixel_data_handlers.apply_voi_lut(file.pixel_array, file) for file in files]
         fig, axes = plt.subplots(nrows=3, ncols=6, figsize=(24,12))
         start = 80
         i = 0
         for i in range(start,start+18):
             img = images[i]
             file = files[i]
             slice_num = i
             # Plot the image
             x = (i-start) // 6
             y = (i-start) \% 6
             for j in range(len(bound_box["slice_number"])):
                 bb_value = bound_box.loc[j]["slice_number"].astype(int)
                 if (bb_value == slice_num):
                         bb_x = bound_box.iloc[j]["x"].astype(float)
                         bb_y = bound_box.iloc[j]["y"].astype(float)
                         width = bound_box.iloc[j]["width"].astype(float)
                         height = bound_box.iloc[j]["height"].astype(float)
                         rect = patches.Rectangle((bb_x, bb_y), width, height, fill=False, edgecolor="blue", linewidth=2)
```


```
axes[x, y].add_patch(rect)
axes[x, y].imshow(img, cmap="bone")
axes[x, y].set_title(f"Slice: {slice_num}", fontsize=14, weight='bold')
axes[x, y].axis('off')
#Each folder is one patient, and the image shown below is in the axial plane
```



The segmentations folder contains segmentations for all slices in a scan in the nifti file format. They were segmented automatically using a 3D UNET model, and radiologists modified and approved the segmentations. However, we need to be careful about the orientation of each.

Please be aware that the NIFTI files consist of segmentation in the sagittal plane, while the DICOM files are in the axial plane.

```
In [12]: #Load segmentation + transpose it to align it with the axial plane image
         seg_path = f"../rsna-2022-cervical-spine-fracture-detection/segmentations/{id_num}.nii"
         ex = nib.load(seg_path)
         seg = ex.get_fdata()
         seg = seg[:, ::-1, ::-1].transpose(2, 1, 0)
         seg.shape
         fig, axes = plt.subplots(nrows=3, ncols=6, figsize=(24,12))
         start = 80
         for i in range(start, start+18):
             mask = seg[i]
             slice_no = i
             # Plot the image
             x = (i-start) // 6
             y = (i-start) \% 6
             axes[x, y].imshow(mask, cmap='inferno')
             axes[x, y].set_title(f"Slice: {slice_no}", fontsize=14, weight='bold')
             axes[x, y].axis('off')
```



In [146... | np.unique(seg[85])

Out[146]: array([0., 3., 4.])

This output of 0,5,6 implies that we are seeing an overlap between C5 and C6

Some notes from research review:

- Utilize a FCN (long and short skip connections speed up the convergence of the learning process) ResNet50 LSTM and HMM
- ResNet50 architecture will feature extract and classify the spine (replace the FC layer of Resnet34 with two FC layers, each being followed by batch normalization and ReLU function)1
 - Initial learning rate 0.01 with mini-batch of 192, reduced in half when no improvement on validation set for 30 epochs
 - The first FC layer reduces 512 to 32 and the second will reduce dimensionality to scalar value
 - Train CNN to make slice-level classification by minimizing binary cross-entropy loss
 - Also will want to use ADAM (0.9, 0.999, 1e-8)
 - Train for 400 epochs
 - Ref
- then an bidirectional LSTM-RNN will run on the resulting probability of the presence of a fracture
 - 256 hidden units use 600 epoch
- Also don't forget to take some of the training for validation